

Reagent Cooling Accessory



**OPTIONAL ACCESSORY FOR ChemWell®
2900 SERIES AUTOMATED ANALYZERS**

**Protects reagents from high ambient
temperatures, ideal for biochemistry,
prolongs the viability of reagents when
in use on-board the analyzer**

- Cooling Capability: 12°C to 15°C below ambient
- Automatic Thermal Control
- Insulated Reagent Racks:
27 position rack standard
44 position rack optional
- Dimensions: 15" (38 cm) L x 6"
(15 cm) W x 8" (20cm) H. Approx.
weight 8 lbs. (3.6 kg)

Reagent Cooling Accessory
Chem Well - RCA (230V) 029040

Reagent Cooling Accessory
Chem Well - RCA (115V) 029041

 **AWARENESS®
TECHNOLOGY, INC.**

sales@awaretech.com

